

# Anti-GCLC Rabbit antibody



Catalog Number: 168667

## Product name

Anti-GCLC Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

GCLC Rabbit monoclonal antibody

## Dilutions

1:1000-1:5000

## Preparation


1:20

Antigen: A synthetic peptide of human GCLC

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of GCLC in C6 cell lysates using GCLC Rabbit mAb(1:1000 diluted).Predicted band size:73kDa.Observed band size:73kDa.

## Clonality

Monoclonal