

Histone H1.2 antibody



Catalog Number: 111367

Product name

Histone H1.2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Histone H1.2 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in Jurkat cells, human testis tissue, L02 cells, MCF7 cells. Observed molecular weight by Western-blot: 32-33 kDa

Preparation

This antibody was obtained by immunization of Histone H1.2 recombinant protein (Accession Number: BC002649). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111367(HIST1H1C antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed MCF-7 cells using Catalog No:111367(Histone H1.2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)