

Anti-Serpina7 antibody



Catalog Number: 102367

Product name

Anti-Serpina7 antibody

Immunogen

[Mouse SerpinA7 \(His Tag\) recombinant protein](#)

Specificity

Mouse Serpin A7 / TBG

Antibody description

Rabbit polyclonal to SerpinA7

Preparation

Produced in rabbits immunized with purified, recombinant Mouse SerpinA7 (rM SerpinA7; NP_808588.3; Met 1-Gln 418). SerpinA7 specific IgG was purified by mouse SerpinA7 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse SERPINA7. The detection limit for Mouse SERPINA7 is approximately 0.078 ng/well.

IHC-P: 0.1-2 µg/mL

Validations



Serpina7 / TBG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse SERPINA7 in mouse liver with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to hepatocyte.