

Human CSNK1G1 (His & GST Tag) recombinant protein



Catalog Number: 502367

General Information

Gene Name Synonym

CKI-gamma 1; EC 2.7.11.1

Protein Construction

A DNA sequence encoding full length of human CSNK1G1 (NP_071331.2) (Met 1-Lys 422) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

The specific activity was determined to be 4 nmol/min/mg using casein as substrate.

Purity

> 94 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human CSNK1G1/GST chimera consists of 659 amino acids and has a calculated molecular mass of 76.3 kDa. It migrates as an approximately 66 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, 10% gly, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CSNK1G1 / CKI-gamma 1 Protein (His & GST Tag) SDS-PAGE