

# PDCD6 antibody



Catalog Number: 113767

## Product name

PDCD6 antibody

## Immunogen

[Human PDCD6 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDCD6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, K-562 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of PDCD6 recombinant protein (Accession Number: NM\_013232). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113767(PDCD6 antibody) at dilution of 1:500



IP Result of anti-PDCD6 (IP:Catalog No:113767, 3ug; Detection:Catalog No:113767 1:500) with HeLa cells lysate 3200ug.