

Anti-phospho-Amyloid Precursor Protein (Thr743) Rabbit antibody



Catalog Number: 167767

Product name

Anti-phospho-Amyloid Precursor Protein (Thr743)
Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

phospho-Amyloid Precursor Protein (Thr743)
Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr743 of human Amyloid Precursor Protein

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of phospho-Amyloid Precursor Protein (Thr743) in HeLa,U2OS,C6 cell lysates using phospho-Amyloid Precursor Protein (Thr743) Rabbit mAb(1:1000 diluted).Predicted band size:87kDa.Observed band size:100kDa.