

Anti-NOLC1 Rabbit antibody



Catalog Number: 168767

Product name

Anti-NOLC1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

NOLC1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NOLC1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000-1/10000

IHC: 1/20

ICC/IF: 1/100

FC: 1/50

Validations



Western blot detection of NOLC1 in Jurkat cell lysates using NOLC1 Rabbit mAb(1:500 diluted).Predicted band size:74kDa.Observed band size:110kDa.