

Mouse Cxcl2/Mip-2/Mip2/Scyb2 recombinant protein



Catalog Number: 518767

General Information

Gene Name Synonym

MIP-2; chemokine ligand 2; C-X-C motif chemokine 2; GRO beta; GRO2; GROB; Gro-beta; Growth-regulated protein beta; Macrophage Inflammatory Protein-2-alpha; melanoma growth stimulatory activity beta; cxcl2; MGSA-b; MGSA-beta; MIP2A; MIP2-alpha; SCYB2.

Protein Construction

A DNA sequence encoding the Mouse Cxcl2/Mip-2/Mip2/Scyb2 Ala28-Asn100 is expressed.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Mouse Cxcl2/Mip-2/Mip2/Scyb2 recombinant protein