

DCXR antibody



Catalog Number: 109767

Product name

DCXR antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DCXR Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human liver tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of DCXR recombinant protein (Accession Number: NM_016286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109767(DCXR antibody) at dilution of 1:500