

NDUFB8 antibody



Catalog Number: 113077

Product name

NDUFB8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDUFB8 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human liver tissue, human heart tissue. Positive WB detected in mouse liver tissue, 3T3-L1 cells, HepG2 cells, human heart tissue, mouse brain tissue, mouse cerebellum tissue, mouse heart tissue, mouse skeletal muscle tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 19-22 kDa

Preparation

This antibody was obtained by immunization of NDUFB8 recombinant protein (Accession Number: NM_005004). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113077(NDUFB8 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:113077(NDUFB8 antibody) at dilution of 1:50 (under 10x lens)


NDUFB8 antibody

Catalog Number: 113077



Immunohistochemical of paraffin-embedded human liver using Catalog No:113077(NDUFB8 antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-NDUFB8 (IP:Catalog No:113077, 4ug; Detection:Catalog No:113077 1:1000) with mouse liver tissue lysate 4000ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:113077(NDUFB8 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

