

# Anti-FBXO6 antibody



Catalog Number: 106077

## Product name

Anti-FBXO6 antibody

## Specificity

Human FBXO6

## Antibody description

Rabbit polyclonal to FBXO6

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human FBXO6, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 10-20 µg/ml

## Validations



Items	Lanes	A	B
Sample (whole cell lysate)		293T	Hela
Sample volume (µg/lane)		30	30
Gel	13%SDS-PAGE reducing gel		
Recommended concentration	10-20µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L) at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 34 kDa Observed band size : 37 kDa		

FBXO6 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot