

# BAX antibody

Catalog Number: 107077

## Product name

BAX antibody

## Specificity

Human; other species not tested.

## Antibody description

BAX Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue, human liver cancer tissue. Positive WB detected in HEK-293 cells, human testis tissue. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of BAX recombinant protein (Accession Number: NM\_138761). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107077(BAX antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107077(BAX Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107077(BAX Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9).

# BAX antibody



Catalog Number: 107077

---

No:107077(BAX Antibody) at dilution of 1:200

(under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).