

# Anti-EGF (5E9) Mouse antibody



Catalog Number: 168077

## Product name

Anti-EGF (5E9) Mouse antibody

## Specificity

Human

## Antibody description

EGF (5E9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant Human EGF protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 20% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1

## Applications


## ELISA

### Validations



sELISA detection of commercial EGF drug by different dilution.

sELISA detection of commercial EGF drug by different dilution.



Standard Curve for Human EGF: Capture Antibody Mouse mAb 168077 (5E9-A10-A6) to Human EGF at 4u03bcg/ml and Detector Antibody Mouse mAb 201165-7C11(7C11-G10-H8) to Human EGF at 1u03bcg/ml.