

Anti-EGF (5E9) Mouse antibody



Catalog Number: 168077

Product name

Anti-EGF (5E9) Mouse antibody

Specificity

Human

Antibody description

EGF (5E9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant Human EGF protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 20% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications


ELISA

Validations



sELISA detection of commercial EGF drug by different dilution.

sELISA detection of commercial EGF drug by different dilution.



Standard Curve for Human EGF: Capture Antibody Mouse mAb 168077 (5E9-A10-A6) to Human EGF at 4u03bcg/ml and Detector Antibody Mouse mAb 201165-7C11(7C11-G10-H8) to Human EGF at 1u03bcg/ml.