

BLOC1S3 Polyclonal Antibody



Catalog Number: 163477

Product name

BLOC1S3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BLOC1S3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-202 of human BLOC1S3 (NP_997715.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - BLOC1S3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using BLOC1S3 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - BLOC1S3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using BLOC1S3 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - BLOC1S3 Polyclonal Antibody

BLOC1S3 Polyclonal Antibody


Catalog Number: 163477

Immunohistochemistry of paraffin-embedded human breast cancer using BLOC1S3 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - BLOC1S3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using BLOC1S3 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - BLOC1S3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse testis using BLOC1S3 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - BLOC1S3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse kidney using BLOC1S3 antibody at dilution



of 1:100 (40x lens).

Immunofluorescence - BLOC1S3 Polyclonal Antibody

Immunofluorescence analysis of NIH-3T3 cells using BLOC1S3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence - BLOC1S3 Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using BLOC1S3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.


BLOC1S3 Polyclonal Antibody

Catalog Number: 163477



Immunofluorescence - BLOC1S3 Polyclonal Antibody


Immunofluorescence analysis of C6 cells using BLOC1S3 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - BLOC1S3 Polyclonal Antibody

Immunofluorescence analysis of NIH-3T3 cells using BLOC1S3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - BLOC1S3 Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using BLOC1S3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.