

SNX22 antibody

Catalog Number: 115477

Product name

SNX22 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SNX22 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of SNX22 recombinant protein (Accession Number: BC019655). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115477(SNX22 Antibody) at dilution of 1:1000



IP Result of anti-SNX22 (IP:Catalog No:115477, 4ug; Detection:Catalog No:115477 1:300) with mouse liver tissue lysate 6400ug.