

IYD antibody

Catalog Number: 111877

Product name

IYD antibody

Specificity

Human; other species not tested.

Antibody description

IYD Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of IYD recombinant protein (Accession Number: NM_203395). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111877(IYD antibody) at dilution of 1:1000