

MSRB3 antibody



Catalog Number: 112877

Product name

MSRB3 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

MSRB3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 26 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of MSRB3 recombinant protein (Accession Number: XM_024448918). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112877(MSRB3 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality