

# MYOT Polyclonal Antibody



Catalog Number: 164877

## Product name

MYOT Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to MYOT

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 259-314 of human MYOT (NP\_001129412.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:10 - 1:100

## Validations



Western blot - MYOT Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using MYOT antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - MYOT Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using MYOT antibody at dilution of



1:100 (40x lens).

# MYOT Polyclonal Antibody



Catalog Number: 164877

---

Immunofluorescence - MYOT Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using MYOT antibody . Blue: DAPI for nuclear staining.