

Adalimumab reference antibody



Catalog Number: 176877

Product name

Adalimumab reference antibody

Specificity

Human TNFSF2 / TNFa

Antibody description

Anti-TNFSF2 / TNFa Antibody (Adalimumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-TNFSF2 / TNFa Antibody (adalimumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFSF2 / TNFa Antibody (adalimumab) is 98.86%, determined by SEC-HPLC.