

DBI antibody



Catalog Number: 109877

Product name

DBI antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DBI Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human colon tissue. Observed molecular weight by Western-blot: 10kd

Preparation

This antibody was obtained by immunization of DBI recombinant protein (Accession Number: NM_020548). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109877(DBI antibody) at dilution of 1:500