

EIF1B antibody

Catalog Number: 110177

Product name

EIF1B antibody

Specificity

Human; other species not tested.

Antibody description

EIF1B Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, HepG2 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 14 kDa

Preparation

This antibody was obtained by immunization of EIF1B recombinant protein (Accession Number: NM_005875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110177(EIF1B antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110177(EIF1B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

EIF1B antibody



Catalog Number: 110177

No:110177(EIF1B Antibody) at dilution of 1:50 (under 40x lens)