

Anti-Fetuin-A antibody



Catalog Number: 106177

Product name

Anti-Fetuin-A antibody

Immunogen

[Mouse Fetuin-A \(His Tag\) recombinant protein](#)

Specificity

Mouse Fetuin-A / AHSG

Antibody description

Rabbit polyclonal to Fetuin-A

Preparation

Produced in rabbits immunized with purified, recombinant Mouse Fetuin-A / AHSG (rM Fetuin-A / AHSG; NP_032037.2; Met1-Ile345). Fetuin-A / AHSG specific IgG was purified by Mouse Fetuin-A / AHSG affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse Fetuin-A / AHSG. The detection limit for Mouse Fetuin-A / AHSG is approximately 0.00245 ng/well.

IHC-P: 0.1-1 µg/mL

Validations



Mouse AHSG Antibody +

Mouse AHSG Antigen -

Fetuin-A / AHSG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse AHSG in mouse liver with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections). The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary antibody and antigen (recombinant protein).