

CRP antibody

Catalog Number: 107177

Product name

CRP antibody

Recommended Dilution:

WB: 1:1000-1:10000


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

CRP Mouse Monoclonal antibody. Positive WB detected in serum from mouse injected with bacteria. Observed molecular weight by Western-blot: 27 kDa




Preparation

This antibody was obtained by immunization of CRP recombinant protein (Accession Number: NM_001329058). Purification method: Protein A purified.

Serum from mouse injected with bacteria were subjected to SDS PAGE followed by western blot with Catalog No:107177 (CRP Antibody) at dilution of 1:2000. Horseradish peroxidase conjugated protein A was used for signal developing instead of secondary antibody.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB

Serum from normal mouse (left) and mouse injected with bacteria (right) were subjected to SDS PAGE followed by western blot with Catalog No:107177 (CRP Antibody) at dilution of 1:500. Horseradish peroxidase conjugated protein A was used for signal developing instead of secondary antibody.

Dilutions