

# Glenzocimab reference antibody



Catalog Number: 177177

## Product name

Glenzocimab reference antibody

## Specificity

Human GP6 / Glycoprotein-6

## Antibody description

Anti-GP6 / Glycoprotein-6 Antibody (Glenzocimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-GP6 / Glycoprotein-6 Antibody (glenzocimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-GP6 / Glycoprotein-6 Antibody (glenzocimab) is 100%, determined by SEC-HPLC.