

# Anti-APG5L/ATG5 (8F1) Mouse antibody



Catalog Number: 168177

## Product name

Anti-APG5L/ATG5 (8F1) Mouse antibody

## Specificity

Human, Rat

## Antibody description

APG5L/ATG5 (8F1) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic peptide.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB


## Dilutions

WB: 1/1000

## Validations



Western blot detection of Atg5 in HCT116, Rat Brain, K562 and C6 cell lysates using Atg5 mouse mAb (1:1000 diluted). Predicted band size: 55KDa. Observed band size: 55KDa.



Western blot detection of Atg5 in Rat Brain and C2C12 cell lysates using Atg5 mouse mAb (1:1000 diluted). Predicted band size: 55KDa. Observed band size: 55KDa.