

PAPOLG antibody

Catalog Number: 113577

Product name

PAPOLG antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

PAPOLG Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in HeLa cells, MCF-7 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 83 kDa

Preparation

This antibody was obtained by immunization of PAPOLG recombinant protein (Accession Number: NM_022894). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113577(PAPOLG Antibody) at dilution of 1:1000



IP Result of anti-PAPOLG (IP:Catalog No:113577, 4ug; Detection:Catalog No:113577 1:1000) with HeLa cells lysate 3200ug.



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog

PAPOLG antibody



Catalog Number: 113577

No:113577(PAPOLG Antibody) at dilution of 1:50 (under 40x lens)