

# MC1R Polyclonal Antibody



Catalog Number: 163977

## Product name

MC1R Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to MC1R

## Preparation

Antigen: Recombinant protein of human MC1R

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:1000 - 1:2000

IHC 1:50 - 1:200

## Validations



## Western blot - MC1R Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MC1R antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.