

# Anti-SIRP alpha/CD172a antibody



Catalog Number: 105977

## Product name

Anti-SIRP alpha/CD172a antibody

## Immunogen

[Mouse SIRP alpha/CD172a \(His Tag\) recombinant protein](#)

## Specificity

Mouse SIRP alpha / CD172a

## Antibody description

Rabbit monoclonal to SIRP alpha/CD172a

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse SIRP alpha / CD172a (rM SIRP alpha / CD172a; BAA20376.1; Met1-Asn373) and conjugated with APC under optimum conditions, the unreacted APC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



SIRP alpha / CD172a Antibody (APC), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse SIRPA(CD172a) expression on BABL/c bone marrow cells. Cells were stained with APC-conjugated anti-Mouse SIRPA(CD172a). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.