

Anti-SIRP alpha/CD172a antibody



Catalog Number: 105977

Product name

Anti-SIRP alpha/CD172a antibody

Immunogen

[Mouse SIRP alpha/CD172a \(His Tag\) recombinant protein](#)

Specificity

Mouse SIRP alpha / CD172a

Antibody description

Rabbit monoclonal to SIRP alpha/CD172a

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse SIRP alpha / CD172a (rM SIRP alpha / CD172a; BAA20376.1; Met1-Asn373) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



SIRP alpha / CD172a Antibody (APC), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse SIRPA(CD172a) expression on BABL/c bone marrow cells. Cells were stained with APC-conjugated anti-Mouse SIRPA(CD172a). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.