

Anti-TLR2 Rabbit antibody



Catalog Number: 166977

Product name

Anti-TLR2 Rabbit antibody

Specificity

Mouse

Antibody description

TLR2 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse TLR2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/20

Validations



Western blot detection of TLR2 in Raw264.7 cell lysates using TLR2 Rabbit pAb(1:1000 diluted). Predicted band size:89KDa. Observed band size:89KDa.