

# Anti-TLR2 Rabbit antibody



Catalog Number: 166977

## Product name

Rabbit IgG

Anti-TLR2 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF

Mouse

## Dilutions

## Antibody description

WB: 1/1000

TLR2 Rabbit polyclonal antibody

IHC: 1/50

## Preparation


ICC/IF: 1/20

Antigen: Recombinant protein of mouse TLR2

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of TLR2 in Raw264.7 cell lysates using TLR2 Rabbit pAb(1:1000 diluted).Predicted band size:89KDa.Observed band size:89KDa.

## Ig Type