

DIS3L2 antibody



Catalog Number: 109977

Product name

DIS3L2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DIS3L2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse embryo tissue. Observed molecular weight by Western-blot: 65 kDa-99 kda

Preparation

This antibody was obtained by immunization of DIS3L2 recombinant protein (Accession Number: BC030113). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109977(DIS3L2 Antibody) at dilution of 1:600