

PRRG2 antibody



Catalog Number: 114277

Product name

PRRG2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PRRG2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of PRRG2 recombinant protein (Accession Number: NM_001316335). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114277(PRRG2 antibody) at dilution of 1:300