

Anti-TAK1 (3G1) Mouse antibody



Catalog Number: 168277

Product name

Anti-TAK1 (3G1) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

TAK1 (3G1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human TAK1 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of TAK1 in K562, HeLa, MCF7, COS7, 3T3 and C6 cell lysates using TAK1 mouse mAb(dilution 1:2000). Predicted band size: 67kDa. Observed band size: 67,78kDa.