

# CHRNA1 antibody



Catalog Number: 109277

## Product name

CHRNA1 antibody

IHC: 1:20-1:200

IF: 1:50-1:500


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

CHRNA1 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Positive IF detected in C2C12 cell. Positive WB detected in mouse heart tissue, rat brain tissue. Observed molecular weight by Western-blot: 54 kDa




## Preparation

This antibody was obtained by immunization of CHRNA1 recombinant protein (Accession Number: BC006314). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109277(CHRNA1 Antibody) at dilution of 1:300

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:109277(CHRNA1 Antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, IHC, WB, IF

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Immunohistochemistry of paraffin-embedded


# CHRNA1 antibody



Catalog Number: 109277

---

human skeletal muscle tissue slide using Catalog No:109277(CHRNA1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed C2C12 cell using Catalog No:109277(CHRNA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)