

IGF1B-Specific antibody



Catalog Number: 111677

Product name

IGF1B-Specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

IGF1B-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human liver cancer.
Positive WB detected in PC-3 cells, human brain tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001111285).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111677(IGF1B-Specific antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:111677(IGF1B-Specific antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:111677(IGF1B-Specific antibody) at dilution of 1:50 (under 40x lens)