

# Anti-BLCAP antibody



Catalog Number: 176677

## Product name

Anti-BLCAP antibody

## Specificity

Human, Rat

## Antibody description

Mouse monoclonal antibody to BLCAP

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 50% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

WB, IHC-P, ICC

## Dilutions

WB: 1:1,000-1:2,000

ICC: 1:50-1:200

IHC-P: 1:50-1:200

## Validations




Fig1: Western blot analysis of BLCAP on recombinant protein tissue lysate using anti-BLCAP antibody at 1/5,000 dilution.




Fig2: ICC staining BLCAP (green) in HepG2 cells. The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig3: Immunohistochemical analysis of paraffin-embedded human breast tissue using anti-BLCAP antibody. Counter stained with hematoxylin.




Fig4: Immunohistochemical analysis of paraffin-embedded rat cervix tissue using anti-BLCAP antibody. Counter stained with hematoxylin.