

C1orf31 antibody



Catalog Number: 108677

Product name

C1orf31 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C1orf31 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, mouse liver tissue. Observed molecular weight by Western-blot: 10-18 kDa

Preparation

This antibody was obtained by immunization of C1orf31 recombinant protein (Accession Number: NM_001012985). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108677(C1orf31 antibody) at dilution of 1:1000



IP Result of anti-C1orf31 (IP:Catalog No:108677, 4ug; Detection:Catalog No:108677 1:500) with HEK-293 cells lysate 1760ug.