

# SLC43A1 antibody



Catalog Number: 115377

## Product name

SLC43A1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SLC43A1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of SLC43A1 recombinant protein (Accession Number: NM\_003627). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: N/A

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115377(SLC43A1 antibody) at dilution of 1:0