

# TTC26 antibody



Catalog Number: 116377

## Product name

TTC26 antibody

## Immunogen

[Human TTC26 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TTC26 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 60-64 kDa

## Preparation

This antibody was obtained by immunization of TTC26 recombinant protein (Accession Number: NM\_024926). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116377(TTC26 Antibody) at dilution of 1:300