

Synaptotagmin-4 antibody



Catalog Number: 115777

Product name

Synaptotagmin-4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Synaptotagmin-4 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, A431 cells, mouse brain tissue. Observed molecular weight by Western-blot: 48-50 kDa

Preparation

This antibody was obtained by immunization of Synaptotagmin-4 recombinant protein (Accession Number: NM_020783). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115777(SYT4 antibody) at dilution of 1:500