

ADAT2 antibody

Catalog Number: 107777

Product name

ADAT2 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

ADAT2 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis tissue. Positive IF detected in A431 cells. Positive WB detected in A431 cells, mouse brain tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 25kd




Preparation

This antibody was obtained by immunization of ADAT2 recombinant protein (Accession Number: BC037955). Purification method: Antigen affinity purified.

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:107777(ADAT2 antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:107777(ADAT2 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF



Dilutions

Recommended Dilution:

WB: 1:200-1:2000


Immunohistochemical of paraffin-embedded

ADAT2 antibody



Catalog Number: 107777

human hepatocirrhosis using Catalog No:107777(ADAT2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A431 cells using Catalog No:107777(ADAT2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)