

Human IL-5/IL5 / Interleukin 5 recombinant protein



Catalog Number: 502087

General Information

Gene Name Synonym

B-cell differentiation factor I; Eosinophil differentiation factor; T-cell replacing factor

Protein Construction

A DNA sequence encoding the human IL5 (P05113) (Ile20-Ser134) was expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Activity

Measured in a cell proliferation assay using TF-1 human erythroleukemic cells.

The ED₅₀ for this effect is typically 0.1-0.5 ng/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human IL5 consists of 115 amino acids and predicts a molecular mass of 13.1 KDa. It migrates as an approximately 27 KDa band in SDS-PAGE under non-reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL5 / Interleukin 5 Protein SDS-PAGE