

# FMO1 Polyclonal Antibody



Catalog Number: 165087

## Product name

FMO1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to FMO1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 363-532 of human FMO1 (NP\_002012.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

## Validations



Western blot - FMO1 Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using FMO1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 150s.



Immunofluorescence - FMO1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using FMO1 antibody. Blue: DAPI for nuclear staining.