

# CAMSAP2 Polyclonal Antibody



Catalog Number: 166087

## Product name

CAMSAP2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CAMSAP2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CAMSAP2 (NP\_982284.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - CAMSAP2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CAMSAP2 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.