

TMED4 antibody



Catalog Number: 116087

Product name

TMED4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMED4 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of TMED4 recombinant protein (Accession Number: NM_182547). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116087(TMED4 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:116087(TMED4 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)