

Anti-Caspase-3 Rabbit antibody



Catalog Number: 167087

Product name

Anti-Caspase-3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Caspase-3 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Caspase-3 p12

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP


Dilutions

WB: 1/500-1/1000


IHC: 1/100-1/500

ICC/IF: 1/50-1/200

Validations




Immunofluorescent analysis of HeLa cells, using



Caspase-3 p12 Antibody.

Western blot analysis of Caspase-3 p12 expression in HeLa cell treated with staurosporine lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using Caspase-3 p12 Antibody.