

Anti-Cytokeratin 5 (7H5) Mouse antibody



Catalog Number: 168087

Product name

Anti-Cytokeratin 5 (7H5) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

Cytokeratin 5 (7H5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human KRT5 beta protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/2000

Validations



Western blot detection of Keratin(50-60kDa) in C6,3T3,A549,MCF7 and HeLa cell lysates using Keratin(50-60kDa) mouse mAb (1:2000 diluted).Predicted band size50~60KDa.Observed band size:50~60KDa.