

# Anti-Biglycan Rabbit antibody



Catalog Number: 169087

## Product name

Anti-Biglycan Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Biglycan Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Biglycan

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Ig Type

Rabbit IgG

## Applications

WB, IP, FC


## Dilutions

WB: 1/500-1/2000

IP: 1/50

FC: 1/100

## Validations



Western blot analysis of Biglycan expression in C2C12 cell lysate.